



Baby Buddy™ Infant CPR Manikin

LF03720U,
LF03721U (5-pk), LF03722U (10-pk)

INSTRUCTION MANUAL



DO NOT STORE IN TEMPERATURES ABOVE 120°

5-Year Warranty

ABOUT THE SIMULATOR

The **Life/form® Baby Buddy™** infant CPR manikin is a realistic and economical manikin for the training of infant CPR skills. Features include lifelike body and head, sanitary airway system, head tilt, realistic chest rise, palpable landmarks, and an easy-to-use, lightweight design.

LIST OF COMPONENTS

- Torso
- Chest
- Head
- Aluminum Screw Post
- 2 Nylon Washers
- Pkg. of 10 Lung Bags
- Lung Installation Tool



SET UP

Your **Baby Buddy™** manikin is shipped disassembled. Simply follow these easy steps to assemble your manikin.



Figure 1



Figure 2

ATTACHING THE HEAD

Position the head over the neck stub of the torso. Insert the screw post through one of the nylon washers and the corresponding holes in the head and neck. (See *Figure 1.*) You may find this easier to do if you turn the manikin over. Securely screw the end onto the screw post over the second washer. (See *Figure 2.*)

Actual product may vary slightly from photo. Nasco reserves the right to change product color, materials, supplies, or function as needed.



Figure 3



Figure 4

INSTALLING THE LUNG BAG

Insert the lung installation tool into the lung bag. (See *Figure 3*.) Fold the round sides of the lung around the tool. Feed the lung into the mouth and through the hole in the back of the head. (See *Figure 4*.)



Figure 5



Figure 6



Figure 7



Figure 8

With the head tilted back, grasp the bag with your other hand and remove the tool. (See *Figure 5*.) Gently pull the rest of the bag through until the barrier flaps reach the mouth. Lay the lung over the corresponding area of the torso. Make sure that the bag is not twisted or pinched, but is open for its entire length. Open the barrier flaps to protect the face. (See *Figure 6*.) Flex the foam chest piece slightly to snap it onto the torso over the lung (See *Figure 7*.) Your **Baby Buddy™** is now ready to use. (See *Figure 8*.)

MAINTENANCE

Used lung bags should be removed and discarded after each use. Always use a new lung bag for each student. The manikin's face may be cleaned with a mild detergent and water if it becomes soiled. Do not attempt to launder the foam pieces or immerse them in water.

SUPPLIES AND REPLACEMENT PARTS

LF03723U **Baby Buddy™** Lung Bags, Pkg. of 100

LF03724U **Baby Buddy™** Carry Bags (Holds 5 **Baby Buddy™** Manikins)

LF03738U **Baby Buddy™** Lung Bag Tools, Pkg. of 4

Other Available **Life/form** Simulators

- LF03623U Infant Airway Management Trainer with Stand
- LF03693U **Basic Buddy**® CPR Manikin
- LF03709U Infant **CRISIS**™ Manikin
- LF06001U CPR Prompt® Adult/Child Manikin
- LF06012U CPR Prompt® Infant Manikin
- LF01420U C.H.A.R.L.I.E.



LF03709U



LF06012U



ONLINE: NascoHealthcare.com
EMAIL: lifeform@enasco.com
PHONE: 1.800.431.4310
ADDRESS: 901 Janesville Ave.
Fort Atkinson, WI 53538